

## Waiting Times Report

Last Modified on 17/04/2024 10:38 am AEST

The Waiting Times Report analyses patient waiting times and can be generated for each doctor and branch in the practice. It's useful for analysing how well doctors work to time and for planning how many doctors need to be scheduled on specific days.

To open the report:

- 1. Go to Zedmed's Reception tab.
- 2. Select Daily Reports > Waiting times.

The Waiting Times Report screen will open.

- 3. Select the version of the report you want to open:
  - **Detailed report** lists every patient seen on the day.
  - Waiting times by hour calculates the average length of time patients waited for each hour of the day.
  - Patients per hour by doctor- calculates the throughput of patients per hour.
- 4. Select the report's Criteria.

You can focus on a specific **Branch** or **Treating Doctor**, and set the date and time range that you want the include.

Use **Waiting Time defined as** to set what actions define the start and end of the period the patient was waiting. See the **How waiting times are measured** section below to learn more.

5. Select **Print** to print out the report, or **Screen** open it. The **Screen** view has print and save menu options.

		o i i i i i i i i i i i i i i i i i i i			of 3 🕨	M	Close			
Branch 1									Page 1	
		١	Nait	ing T	imes R	eport				
Criteria						1,541	Group By			
Date Range 29/07/2020 - 2/08/2021 Branch Branch 1						None				
	Treating Doctor All Waiting Time Appt/Attend (later of) to Admission time									
1000		20 20				- F	Order By Attendance Time			
Date	Patient	Brn	Doc	Appt	Attended	Admittee	d Billed	Waited		
5/07/2021	Holaway, Knut	MED	PD	08:30	08:47	08:47	08:54	00:00		
1/07/2021	Star, Mike	MED	PD	09:00	09:03	10:32	10:33	01:29		
9/07/2021	Holaway, Knut	MED	PD	00:00	09:12	09:12	09:26	00:00		
9/07/2021	Holaway, Knut	MED	PD	00:00	09:12	09:12	09:28	00:00		
9/07/2021	Holaway, Knut	MED	PD	00:00	09:12	11:14	14:09	02:02		
8/07/2021	Holaway, Knut	Waiting Times Report								
8/07/2021	Allison 1, Anna	Report Criteria								
8/07/2021	Holaway, Knut	Detailed rep					aa /az /aa aa			
8/07/2021	Bedker, Irene					Entry	Date	29/07/2020 🚔 📘	02/08/2021	¥
1/07/2021	Holaway, Knut	Group By	None		~	Brand	h	Branch 1		
30/06/202	Barker, Emily	Order By	Atten	dance Time	• •	Treat	ing Doctor	All		
5/07/2021	Barker, Emily		12							
7/06/2021	Black, Emma	O Waiting times by hour				Start	Time	8:00 AM V End	d Time 6:00 Pt	M N
7/06/2021	Black, Emma	O Patient per hour by doctor								
7/06/2021	Black, Emma					Waitin	Waiting Time defined as			
29/07/202	Holaway, Knut					A	ppt/Attend (lat	erof) v to A	Admission time	
19 M	Holaway, Knut									

## How waiting times are measured

The waiting time is controlled by the **Waiting Time defined as** option and will reflect whether the patient is **Admitted**, which populates when the doctor opens the patient's record from the Clinical Waiting Room. If there is no **Admitted** time then the billing time will be used to calculate the patient's exit time instead.

The start of **Waiting Time** can be defined in terms of the booked **Appointment Time** or the **Arrival Time**, normally by the latter of the two values, so a patient that arrives very early for their booked appointment does not compromise the waiting time average for that hour.

Attendance records created as a result of **Billing** via **Patient Details**, which create attendances records for 6.00 am, are excluded from the **Waiting Time** analysis reports as they would be misleading.