

Spellcheck

Last Modified on 21/02/2024 12:42 pm AEDT

Use Zedmed's spell checker to review and check your notes from the Tools menu or the Current Encounter menu. By default, all text you enter will be checked for spelling errors, however, opening the Spell Check allows you to review and action all recommendations.

To run a Spell Check from the encounter menu:

1. Go to Zedmed's **Reception** tab.
2. Select **Clinical Records**.


The **Zedmed Clinical** console will open.

3. Start an encounter with the patient.

For a detailed explanation of how to do this, see the [Start an Encounter](#) article.

When the encounter starts, the **Current Encounter** pane will open and display the clinical modules.

4. Select **Spell Check** from the **Current Encounter** menu.



The Spelling dialog will open and check any text you have entered.

5. Use the Spelling dialog option to ignore or change the spelling.

