

# Summary Views Medications

Last Modified on 22/07/2024 9:51 am AEST


Use the **Medications** tab in **Summary Views** to review a patient's past and current prescriptions.

- Expanding a medication will show the prescriptions.
- Selecting the **Include MHR** icon will show prescriptions uploaded to MHR by other providers.
- Red prescription text means a renewal is due.

Right-clicking a medication or prescription displays additional options including:

- **Filter History View to these Encounter Segments** - shows only and all instances of this medication being prescribed.
- **Re-Prescribe Drug** - opens the dosage screen with the same medication already selected.
- **Prescribe New Strength** - searches MIMS for the medication name in order to allow prescribing of a new strength.
- **Make Non-Current** - Updates the drug's status to one that the patient is not currently taking.
- **Edit drug details** - Opens the **Script Details** screen. Changes to the script will be used when re-prescribing.
- **Delete this drug** - Deletes the drug and all prescriptions for it. This will not delete a drug with electronic prescriptions.

See the [Reprints and Cancellations](#) guide to learn more about medication management.



The screenshot displays the 'Summary Views' interface for medications. On the left is a navigation menu with options like Problems, Referrals, Results, Immunisations, Allergies, Images/ECG, Pregnancies, Documents, Attachments, Measurements, Medications, and Incoming Documents. The main area shows a list of medications with filters for 'Prescribed >=' (set to 'All') and 'Include MHR'. The 'Sort' dropdown is set to 'Medication Text'. The list includes:

- Current Medications**
  - Regular Use
    - Olanzapine 10mg Coated Tablet 10mg [28] - TAKE 1 Q4H TAB AFTER MEALS BID AS DIRECTED
    - Paracetamol 500mg - Ibuprofen 200mg Oral Capsule (500 mg/200 mg) [201] - TAKE 1 BID
    - 23/02/2024 TAKE 1 BID
  - PRN
    - Short Term
- Non-current Medications**
  - Temazepam 10mg Oral Tablet 10mg [25] - TAK
  - 23/02/2024 TAKE 41 NOCTE

A right-click context menu is open over the Paracetamol/Ibuprofen medication, showing the following options:

- Filter History View to these Encounter Segments
- Re-Prescribe Drug...
- Prescribe New Strength...
- Add New
- Edit Drug Details...
- Make Non-Current
- Hide/Archive if not current
- Search based on this drug
- Make All Non-Current
- Make All Current
- Delete this Drug...
- Show Notifications
- Show drugs

Below the medication list is a 'Reference View' section showing 'Allergies: Penicillins' and 'Warnings: ...'.