

Add problem from ICPC

Version: 1.00 | Last Modified on 29/04/2025 8:26 am AEST

Use the Problems module to record an encounter's diagnosis using the ICPC list. ICPC is an International Classification method for Primary Care encounters.

Alternatively, you can record a problem by manually entering it or selecting one from the existing or recent problems list. The method used is determined by the doctor or the practice's guidelines.

Requirements

ICPC requires an active subscription to ICPC-2 Plus.

To add a problem from the ICPC list:

1. Start an encounter with the patient.

For a detailed explanation of how to do this, see the Start an encounter article.

When the encounter starts, the Current Encounter section will open and display the clinical modules.

2. Select **Problems** from the **Current Encounter** menu.

Curr	Current Encounter - Duration: 00:03:01													
	Dems	⊖⊗ _{Drugs}	Observations	Referrals	پنې Immunisations	Allergies	A Diagrams	Templates	CDM		Capture	×		
Undefined Problem														
Θ	RFE:	OC - Cus	tom RFE									Use Coded P	RFEs Manageme	ent Plan
\sim	Cons	ultation	notes											

The Define Problem/REF screen will open (screenshot below).

- 3. In the Problems section, select Predictive Search (ICPS and Common) or ICPC Search.
 - For Predictive Search, enter the problem's name in the Filter Current List field to display the problems to choose from.
 - For ICPC Search, use the Keyword, Code and Chapter fields provided to display the problems to choose from.
- 4. Select the problem.
- 5. As an optional step, you can use the **Extra Text** field to add a note, for example, "left side" or "Head".
- 6. Select the **Select Problem** button (You can also double-click the problem).

This adds the problem to the tab above RFE that previously said Undefined.

To add additional problems, select the blue + icon on the top left and repeat the steps above.

Optional selections:

- Add to History adds the problem to the patient's history.
- Add to Common adds the problem to the practice's Common Problems list.

• **Onset Date** records when the problem started.

This is required for My Health Record Shared Health Summary uploads.

- Accurate To notes the approximation of the onset date.
- Criticality can be used to sort a patient's problems in Summary Views.
- 7. Select **Close** to save the problem.

When you select Refresh, the problem will appear in the patient's **Summary Views > Problems** tab.

Current Encounter						
Head Cold						
RFE:					Use Coded RFEs	Management Plan
						^ ~
		Head Cold (Exis	sting Problem)			
Problem/RFE Selection						
Problems						
Selection Method			_			
O Existing Problems O My	Recent Problems Predictive S	Search (ICPC and Common) O ICP	C Search O Common Problems	O Free Text	Problems / Syn	mptoms Only
Filter						
Filter Current List head						
blackheads cervicogenic headache	head injury head injury with fracture	headache heaviness of the head	symptom or complaint of the heat tension headache	be		
dosed head injury	head injury without fracture	lightheadedness	vascular headache			
duster headache	head lice infestation	migraine headache	whiteheads			
congestion in head head cold	head pain head spinning	post traumatic headache sinus headache				
nedu colu	nead spinning					
Problem	Extra Text					
Head cold		ICPC Code: R74 004		Select Problem	Add to History	Add to Common
✓ Include in Correspondence Problem Details for - Hes	ad Cold (Existing Problem)					
Onset Date 05/10/2021	🗾 Accurate To: 🖲 Day 🔿 Mo	nth OYear O0ther	Criticality Important	~	Left Bight	Current
Comments and Management PI	lan					
						~
						~
				? Help	i d	llose