

Surgical Partners

Last Modified on 22/05/2024 3:12 pm AEST

Surgical Partners is a financial management platform that integrates Zedmed with Xero financial software. This allows you to use the financial information in Zedmed for a range of financial applications from debtors' management to banking and payroll. This functionality is enabled by registering with Surgical Partners and then configuring Zedmed to send Surgical Partners the required data each night.

This is a paid service managed by Surgical Partners. To learn more, see https://surgicalpartners.com.au.

When Zedmed makes its initial connection to Surgical Partners, it will send unpaid and part-paid invoice information first to establish the opening balances for each payer. It will then do a full synchronisation at the scheduled time, which is 4 am by default.

Configuration requirements:

- A Surgical Partners token.
- A user with administrator access to your Zedmed server.

To configure Zedmed's Surgical Partners connection:

- 1. Login into Zedmed with administrator access.
- 2. Select the Management tab.
- 3. Select Practice Setup > Practice.

The Practice Details screen will open.

- 4. Select the Integrations tab.
- 5. In the Surgical Partners section, select the Connect button.

The Surgical Partners Connection dialog will open.

gement Reports Tools Pr	actice Setup WD Setup I Hilities	Branch 1		<u> </u>	
Reception	Management	Reports	Utilities		
Claims		Practice Details			
Debtors	Practice Surgica Branches This pu	l Partners ractice is not connected to Surgical Partners	3 Connect		
3rd Party Batc 2	Bank Accounts Departments Integrations Electro	nchronise data with Surgical Partners	Daily upload time (hour) 4 ical Partners Connection	Manual Synchronisation	
Reports	Claiming Configuration	able Electronic Document Transfer	· []		
Practice Setup			Connect Don't Connect		
WP Setup					
Messages					
AIR					
Stock Management					
Message o edmed Support h: 1300 933 000 :: support@zedmed.com.au V: www.zedmed.com.au	a da anti-				
				Close Cance	

6. Copy the token provided by Surgical Partners into the **Token** field.

Caution: Before selecting the Connect button, read the remaining steps and the Considerations below.

7. Select Connect.

A dialog will ask if you want to start the synchronisation with Surgical Partners.

- If you select Yes Zedmed will send unpaid and part-paid invoice information to Surgical Partners.
- If you select **No** the connection will be established but no data will be transferred.

Confirm		×
?	Do you want to start the synchronisation with Surgical F	Partners now?

8. Review the settings and select Close when you have finished.

Considerations

An initial upload should be performed by selecting **Yes** to the dialog above to upload the unpaid and part-paid invoices. This will allow you to test the connection and create opening balances for the payers. You can then leave the scheduled upload to perform the full synchronization.

You can use the Manual Synchronisation button to force a synchronisation at any time.

To stop the scheduled synchronisation from running, untick **Synchronise data with Surgical Partners**.