

ReturnToWork SA

Last Modified on 22/04/2024 3:15 pm AEST

This process has been replaced by https://help.zedmed.com.au/help/rtwsa-healthlink-smartform

Zedmed supports the ReturnToWork SA Electronic Work Capacity Certificate (eWCC), which makes it easy for doctors to certify patient capacity and prescribe treatment. Zedmed streamlines the process of creating and sending an eWCC, and once you have created a certificate, it can be cloned by other doctors at the same practice to create new certificates for the selected patient.

To create a certificate, see our Clinical eWCC guide.

To enable the eWCC feature:

- 1. Go to Zedmed's **Reception** tab.
- 2. Select Clinical Records.

The Clinical Records screen will open.

3. Select Tools > Global Options.

The Global Options screen will open on the General tab.

4. Select the Enable Return to Work SA electronic Work Capacity Certificate tick box.

You can now use an eWCC to certify patient capacity and prescribe treatment.

| eneral | Communications | Letters & Tem | plates Scanning | Results Drugs | ACIR Immunis | tion Defaults | Chronic Disease Manager | ient Library | |
|---|--|---|--|--|------------------------|--|---|---|--|
| Site Type | | General Pra | ctice | ~ | | | | | |
| Dis Alk Usi Alk Alk Alk Clinica | pensary Active ow Pre/Post Dated S cure Template Edito e ICPC Natural Lang ow Free Text Proble FE or Problem Name I Record Privacy Sel able user group rest | Scripts r Juage ms Required tting riction to view p | Show Medicare I Use Clinical Billin Use Chronic Disu Use Coded RFEs Allow Fore Text Allow Forward D | Expiry on Script g Items Wizard ease Module Problems in 'My I ates for Referra | Recent Problems' Is | Warn or * Notes on Encounter Immunisat Summary N | n Failed HL7 referrals RFE or Problem Name Requ tabs created as a result of jons or by re-prescribing m Views are exempt. | ired recording Allergies, edications from the | |
| Defau Da | lt Problem Onset Ac y C | curate To) Month | () Year | OF | ree Text | | | | |
| Pa | tient Education Link Bulk Bill Show Path. Details Show Attachments Iministrator Override | W:\VM_Setu Form Tab | up_Stuff\Patient lea | ftlet\Murtagh\pe | -index.htm 😰 | | | | |
| Use On Offline | e Practice Postcode Iy use eClinic for pa | for Childhood I thology | mmunisation Sched | ıle | | | | | |
| | RS Server Folder able Return to Work | \\DVICQVID | EV2019B\ZedmedSe Nork Capacity Certi | rver\ ficate | * | | | | |
| | | | | | | | | | |
| | | | | | | | | | |