

# June 25th eRx requirements

Last Modified on 03/06/2024 9:56 am AEST

The `erx.com.au :8830337559` certificate is [expiring on 25 June 2024](#) and the eRx Adapter must be updated to ensure continued operation of the eRx service.

HTTPS will use a different port and the eRx web address will change to HTTPS. If you use a firewall with port/internet restrictions, you may need to open a new port to support the new endpoint.

**Zedmed v36.6.8 will update your eRx certificate and change the communication endpoint from HTTP to HTTPS.**

Learn more about [Zedmed v36.6.8](#)

---