

Zedmed v38.1.2

Version: 1.00 | Last Modified on 08/04/2025 4:03 pm AEST

Clinical	ZED-10187	Fixed error: Failed to create Edge browser control This error occurred when opening Zedmed Clinical on a Terminal Server, which prevented Heidi from being available if multiple users were running Clinical on the same computer.
----------	-----------	--


Zedmed v38.1.1

Clinical	ZED-10164	AV when selecting Clinical Details in Office and Heidi is not enabled Fixed an AV error that can occur in some instances when the patient history is refreshed.
Office	ZED-10160	Create an installer that installs WebView2 automatically Each time Office starts up it will check if WebView2 is installed, if not it will run the installer silently. Logging added at Debug level in the Patients.log which shows the check and if the installer was run.

Zedmed v38.1.0

NEW Prevent specific patients from booking in OLA

Practices can now block existing patients from booking appointments using OLA. The feature is enabled using a tick box in the Patient Details and applies to existing patients with OLA accounts that have been reconciled to their patient record.



The screenshot shows the 'Patient Details' form for Chris Armstrong. The form is divided into several sections: Personal, Practice, Family, Payers, Referrals, Other Contacts, and Alerts. The 'Personal' section includes fields for Family Name, Given Name, Preferred Name, Title, Sex at Birth, Pronouns, Date of Birth, Age, Home address, Suburb, Preferred Contact Method, Mobile phone number, Email, Usual Dr, Medicare, Veteran, Health Care, Safety Net, Pen. Stat., ATSI, and Ethnicity. The 'Practice' section includes File#, Hosp. UR No, Usual Branch, First In, Last In, Type, Status, and Clinical Warnings. The 'Family' section includes Family Name and Family Links. The 'Payers' section includes Fund and No fund set for this patient. The 'Referrals' section includes Referrals and Add/Details buttons. The 'Other Contacts' section includes NOK Name, Relationship, Mobile, and Emerg. Contact Name. The 'Alerts' section includes Alerts. The 'Block patient from online bookings' checkbox is highlighted with a red box.

This is what the patient will see in OLA,

Confirm My Booking

There was a problem booking your appointment, please contact the centre

Patient Selection

Please select the patient attending this appointment from the list below, or select "Lookup Patient" and fill in their details.


☐ Matt

☐ Matilda

☒ Chris

☐ Lookup Patient

Continue



Dr Phillip Davis
FRACGP, MBBS, BBiomedSci (Hons), Clin Dip
Pall Med

Appointment Details

CENTRE:

Medical One - QV

TYPE:

Std Consultation

Select New Appointment

Update

Clinical	ZED-10142	Cannot focus on a disabled or invisible windows Fixed an issue where an error was displayed when opening a patient from the clinical waiting room.
----------	-----------	---

Zedmed v38.0.0

NEW Heidi AI

Originally released as pilot build v37.9.0

Zedmed's Heidi AI integration uses voice input to generate AI-enhanced clinical notes and consultation transcripts. Before it is added to a patient's file, all generated content can be reviewed and manually edited.

To use Heidi AI, each practitioner needs their own paid subscription. To subscribe, select the **Set up Heidi** link in Clinical's Heidi AI widget. You can also use Heidi's 30-day free trial via this link.

To learn more, see the [Heidi AI guide](#).

Scan

MEMOS

Address book

Travel

Batch letter

To do list

ZEDMED | CLINICAL

Current Encounter - Duration: 00:02:04

Undefined Problem

RFE

Subjective: Dizziness, nausea, lightheadedness

Objective: BP: 122/67
Pulse: 62

Assessment & Plan: 1. Dizziness, nausea, lightheadedness
- Assessment: Vertigo
- Investigations: Consider further workup if symptoms persist
- Treatment: Symptomatic management, advise adequate hydration
- Follow-up: Review if symptoms worsen or persist

History View

1 2 3

Search History

Billing Wizard

Oldest Encounter: 03/01/2025

Visible: Last 10 Consultations

03/01/2025 (Fri) 10:27am with Dr Phillip Davis at MED, for 1m 22s

Subjective: Dizziness, nausea, lightheadedness

Objective: BP: 122/67
Pulse: 62

Assessment & Plan: 1. Dizziness, nausea, lightheadedness
- Assessment: Vertigo
- Investigations: Consider further workup if symptoms persist
- Treatment: Symptomatic management, advise adequate hydration
- Follow-up: Review if symptoms worsen or persist

Context

Transcript

Note

General Pr...

Dr Goldlocks

Subjective: Dizziness, nausea, lig

Brief

Detailed

Goldlocks

Super Detailed

Objective: BP: 122/67
Pulse: 62

Assessment & Plan: 1. Dizziness, nausea, lightheadedness
- Assessment: Vertigo
- Investigations: Consider further workup if symptoms persist
- Treatment: Symptomatic management, advise adequate hydration
- Follow-up: Review if symptoms worsen or persist

Review your note before use to ensure it accurately represents the visit.

Push note

NEW OLA document requests

Practices can now ask patients to upload documents (such as referrals) when making online appointments. The document requirement are enabled for specific Appointment Types, and the documents are delivered to the practice mailbox.

To learn more, see the [OLA document requests guide](#).

Screenshot of a Patient being prompted to provide a referral in OLA.

Patient Selection
Please select the patient attending this appointment from the list below, or select "Lookup Patient" and fill in their details.

☒ Matt
☐ Lookup Patient

Please upload any relevant documentation (referral letter, test result etc.) or select "Skip" if you would prefer to bring a copy with you to the appointment.

Drag & drop your referral letter here or click to select your file

[Skip](#) [Continue](#)

Appointment Details
CENTRE: [Redacted]
Your local time (AEDT)
DATE: Wednesday 19 February 2025
TIME: 12:00 pm
TYPE: Std Consultation
[Cancel Appointment](#)

NEW Schedule conflicts preview

When adding or editing a schedule, you can now view conflicts on a dedicated pane, making it easier to review all conflicts and their details in one place. Conflicts often occur when changes are made, and those conflicts will be displayed on the right pane with a red border.

To learn more, see the [Create a Schedule guide](#).

This new feature will allow you to review and modify conflicts

Set Schedules for Doctor or Resource

Doctor: Davis, Phillip Resource: [Dropdown]

Schedules

Start Date	End Date	Cycle (Days)	Branch
1/01/2017		7	Branch 1
1/08/2016		7	Apostro'sand
1/08/2016		7	Branch 3
1/08/2016		7	Branch 4

[Save](#) [Add](#) [Copy](#) [Delete](#) [Graphic Display](#) [Show all schedules](#) [Clear...](#)

Sessions

Day #	Day Of Week	Start Time	End Time	Slot Length	Activity Type	Available Online
2	Tuesday	09:30	12:00	15 Min	Consultation	<input type="checkbox"/>
1	Monday	08:00	12:00	15	Consultation	
2	Tuesday	09:30	12:00	15	Consultation	
2	Tuesday	13:00	16:00	15	Consultation	
3	Wednesday	08:00	12:00	15	Consultation	
4	Thursday	08:00	12:00	15	Consultation	
5	Friday	08:00	12:00	15	Consultation	

[Save](#) [Messages](#) [Add](#) [Delete](#) [Duplicate](#) [Check For Conflicts](#)

Please resolve the following conflicts to save and close the schedule.

A session from Monday 14:00 to 16:00 for a Schedule with Start Date 1/08/2016 for Doctor Phillip Davis conflicts with a session from Monday 13:00 to 15:00 for a Schedule with Start Date 1/01/2017

A session from Tuesday 14:00 to 16:00 for a Schedule with Start Date 1/08/2016 for Doctor Phillip Davis conflicts with a session from Tuesday 13:00 to 16:00 for a Schedule with Start Date 1/01/2017

This practitioner works at 4 branches and has overlaps between 2 branches on Monday and Tuesday. A conflicting session will need to be modified or removed for each overlap.

[Close](#) [Cancel](#)


NEW Tyro Terminal selector

Originally released as pilot build v37.8.2



Currently, for Zedmed Cloud, users are linked to Tyro terminals, and each user can only be linked to one. This causes problems if users move desks or need to use a different terminal.




This new feature allows users to select the required terminal when they log in. For on-premise customers, workstations are linked to Tyro terminals (not users) so no change is needed.





To learn more, see the [Tyro Terminal selector guide](#).




Enhancements

OLA		<p>Practices can set the order of doctors in OLA.</p> <p>The order in which doctors appear in OLA can be modified in Branches > Online Appointment Details. If this feature is used, it overrides the current order, which is determined by the left-to-right doctor order in the Appointment book. To learn more, see the OLA setup guide.</p> 
Office	ZED-9831	<p>A patient's next appointment is now shown on their patient record.</p> <p>Appointments in the past are not shown and will display "No upcoming appointments" if the patient has no future appointments.</p> 

Office	ZED-9846	<p>Unassigned Tasks will no longer show in Task Manager.</p> <p>Task Manager now displays tasks that have not been assigned to staff when selecting the ALL Assignee option. NOT ASSIGNED has been added to the list of Assignees to enable viewing just unassigned tasks.</p>
Office	ZED-10000	<p>In Patient details, the Notepad timestamp now has both the date and time</p> <p>Added time to the timestamp note in patient details, which previously only had the date.</p> 
Clinical	ZED-9633	<p>Secure Send can now be used in Batch Letter Writer for individual letters.</p> <p>Because a batch is sent at once and there is nowhere to add individual emails, Secure Send is provided as a right-click option for each letter.</p> 
Clinical	ZED-9745	<p>Summary Views now supports more file types.</p> <p>Image files can now be attached to the ECG/Images, Incoming Documents and Results tabs. The supported file list is now: jfif, pjpeg, jpeg, jpg, png, pdf.</p> 
Clinical	ZED-9901	<p>Summary Views now support the Docx file type.</p> <p>The docx files can now be attached to Incoming Documents, Documents, Results, and ECG/Images.</p>

Clinical	ZED-9808	<p>'Results' in Summary Views can be sorted by Date Collected.</p> <p>A new option to sort by 'Date Collected' has been added to the Results tab on the Summary Views. When sorting by this option, the appropriate date is displayed for each result listed.</p> 
Clinical	ZED-9872	<p>Results Inbox 'Common Phrases' now support more characters.</p> <p>The length Common Phrases used to update a result in the Results Inbox in Clinical Records has been extended to 320 characters.</p> 
Clinical	ZED-9893	<p>The Script Report now shows all drugs for a patient, even if they changed their name.</p> <p>The Clinical Script Report now filters patients using their unique patient IDs. This update resolves the issue where patients who have changed their names and received medication under different names were not correctly identified.</p> 
Clinical	ZED-9930	<p>Moved the 'Send to Patient button' in the Results Inbox.</p> <p>The Send To Patient button has been moved to the left of the configurable buttons in the Results Inbox. It was previously to the right of the Show Result Complete.</p> 

Office		<p>Updated the Schedule Details Graphic display:</p> <ol style="list-style-type: none"> 1. Removed the Doctor and Resource radio buttons. 2. Changed the Doctor/Resource dropdown to be a label that displays the selected schedule 3. The End Date check box defaults to the current date rather than the scheduled start date. 
--------	--	---

Fixes

Office	ZED-6353	<p>Bulk patient Merge issues</p> <p>The following minor usability issues have been fixed on the Patient Merger - Bulk merge</p> <ol style="list-style-type: none"> 1. Do not proceed with merge if the user clicks on the cross icon on the confirmation dialog. 2. Do not show an error message when the user cancels printing.
Office	ZED-9770	<p>3rd party Batch Report - The description overlapped the GST</p> <p>Long item descriptions no longer obscure the GST Amount in the 3rd Party Batch Report.</p>
Office	ZED-9845	<p>Appointment Reminders - Automatic Reminder SMS is not listing resources</p> <p>When sending automatic reminders and confirmations, the Resource name is now displayed if the appointment is with the resource and not a doctor.</p>
Office	ZED-9847	<p>Quote Not Showing GAP Fee</p> <p>Fixed an issue with Billing and Quoting where the Gap was not included on the invoice/quotation even though the print gap option was selected.</p>
Office	ZED-9867	<p>Account Payer Allocation - Unable to view Invoice error " Print Duplicate error"</p> <p>Third Party Account Enquiry/Allocation screen has been fixed to prevent an error when viewing third party invoices.</p>
Office	ZED-9869	<p>Merge Patient files - MyMedicare registration Fields Are Not Saved After Merging</p> <p>When merging patients with My Medicare details, Zedmed now allows the user to select how these are saved in the merged record.</p>

Office	ZED-9891	Appointment grid "Could not convert variant of type (Null) into type (Integer)" error Display a more user friendly message when trying to attend and appointment has been deleted by another staff member. Message is "The attendance could not be saved, it has been deleted by another member of staff".
Office	ZED-9834	Assertion failure error for third-party payers Office now blocks attempts to use Secure Send for Third-Party Account Payers, preventing errors that previously occurred when this unsupported feature was incorrectly available.
Office	ZED-9947	Unable to print the Referral Doctor label when using the internal Mail Merge template, getting a Print Referral Doctor Label Error Fixed the issue where an error is reported when there are no referrals relating to the Referral Doctor, using the internal mail merge templates.
Office	ZED-9966	Next Available Search function is not updated after clearing a schedule conflict We have resolved an issue in Zedmed Office where appointment slots were being lost due to scheduling conflicts in doctor schedules.
Office	ZED-10009	Update to Waiting Room auto-refresh frequency branch option Values entered that are 15 or less will now be a refresh in minutes time. Values above 15 will be a refresh time in seconds. The label will automatically change from "(seconds)" to "(minutes)". Label caption at top left of waiting room window will also change.
Office	ZED-10028	Cannot send ZedSMS from new branch - Error - 103: Invalid Config Values When creating a new branch, Zedmed will now correctly save the SMS service provider settings. Previously, users could not send SMS messages from newly created branches.
Office	ZED-10071	Remove the Tip of the Day Editor menu option in Office The Tip of the Day feature is no longer supported and has been removed from Office.
Clinical	ZED-8763	To-Do List freeze when open patient file from result inbox Fixed the issue when clicking Action Performed on a To Do List item for a patient opened from the Results Inbox, Clinical appears to hang.
Clinical	ZED-9629	Results Inbox Access Violation when using arrow keys to navigate from a result that contains an image Fixed an error that can occur when using the keyboard to scroll through results in the results inbox.
Clinical	ZED-9676	AV error when opening the results inbox after importing results Fixed an error that can occur when importing results with a patient record open then opening the results inbox.
Clinical	ZED-9732	Rename option on Pathology/Radiology referrals does not work Removed non-functional rename option from the right click menu of pathology and radiology referrals in the Encounter History view.

Clinical	ZED-9733	AV error when deleting a referral using the contextual menu in History View Disabled non-functional delete option in right click menu for referrals shown in the encounter history. Use the summary view section to manage these documents.
Clinical	ZED-9767	Clinical encounter - Duplicate addendum notes Fix issue that was causing addendums to clinical notes to appear duplicated if a second addendum was added.
Clinical	ZED-9768	Access Violation when viewing/opening images from history view. Clinical no longer throws an error when right clicking on an image in the encounter note.
Clinical	ZED-9769	Referral Templates via Encounter - referral module the global template appears as default instead of the My Option default template. The Clinical Referrals dialog has been fixed to ensure we retain the correct template selections.
Clinical	ZED-9771	Referrals made through HealthLink Smartforms still appear in the list of incomplete referrals even after they have been marked as closed. The Outstanding and Requested referral reports in Clinical have been updated to work with HealthLink referrals.
Clinical	ZED-9804	Update Pregnancy Calculator to current year and beyond. The pregnancy calculator has been updated to allow start dates in the range of 2023 to 2043.
Clinical	ZED-9815	Error -- No Measurement Type exists for measurement SPO2. Fixed error "No Measurement Type exists for measurement SPO2" that occurs when trying to use the SPO2 measurement in clinical observations for some users.
Clinical	ZED-9816	Issue with Distributing Letter from Summary View. Clinical letter writer has been fixed where selecting the distribute button would not show the Distribution screen unless the letter had been modified first.
Clinical	ZED-9868	Range Check Error - Result Inbox. Fixed issue that would cause a range check error to be displayed when viewing some results and selecting the Enhanced View tab.
Clinical	ZED-9892	Billing Wizard - Billing details The "Send Medicare Bulk Billing Consent" button in the billing wizard is covering up text in the payer details. The Medicare Bulk Billing Consent button on the Clinical billing wizard have been repositioned to ensure that the patient's fund details are not cut off.
Clinical	ZED-9912	Clinical Results Inbox filters - on Document or Result does not work after a refresh. The Type filter in Results Inbox was not taking effect on a Refresh. This issue has been resolved.

Clinical	ZED-9931	Printing in progress error. Fixed an issue with printing results, where clicking the print button too many times results in an error.
Clinical	ZED-9968	Referral letter writer screen - Printer in Progress error which leads to an AV. The Clinical letter writer will now handle accidental invocations of the printer when it is already busy printing. Previously, this would trigger an error message.
Clinical	ZED-9977	Access Violation when printing two copies of a referral letter. Clinical Letter Writer has been enhanced to prevent an error when a user closes the window while it is printing two copies of the letter.
Clinical	ZED-9993	Creating a new problem using the same code as an existing or resolved problem incorrectly updates the existing record. It is now possible to create new patient problem(s) for a problem code that has already been created.
OLA	ZED-10044	<i>Released in Zedmed v37.8.3 pilot</i> Resolved performance issues related to the Next Available Appointment query <ul style="list-style-type: none"> • Added a configuration to allow the Next available appointment slot request to cache the next available appointment information. This speeds up the response time. • Added a configuration to modify the number of days OLA will look forward for next available appointments. This configuration is managed by Zedmed support and cannot be done in the Zedmed UI.

See all [Zedmed v37 release notes](#)

See all [Zedmed v36 release notes](#)