

# Telehealth on-prem

Version: 1.00 | Last Modified on 08/08/2025 3:45 pm AEST

Conduct remote video consultations with Zedmed patients using an integrated workflow of appointment types and notifications.

Click here to view a one-page [quickstart guide](#)

## Workflow overview

The default workflow ensures the patient joins the Telehealth session before the practitioner. The patient is **Admitted** to the Waiting Room when they tap Join on the Telehealth page on their phone, which starts their video connection. When the practitioner **Attends** the patient from the Waiting Room, the Current Encounter starts and displays the menu the practitioner uses to join Telehealth themselves. If the practitioner opens the patient record manually, a red banner will advise that the patient has joined, but reception will not see an attend time for this consultation.

### Workflow summary:

1. The patient books a Telehealth appointment and receives a confirmation SMS.
2. Before the appointment (e.g., 30 minutes), the patient receives an SMS with the Telehealth link.
3. The patient taps the link to open Telehealth, then taps **Join** (selectable 30 minutes before appointment).
  - The video connection opens on the patient's phone. Displays - 'Your practitioner will arrive shortly'.
  - The patient is automatically **Attended** to the **Waiting Room**.
4. The practitioner sees the patient in the **Waiting Room** and selects **Admit**.

This opens the patient's clinical record and starts the Current Encounter in Zedmed.
5. The practitioner selects Telehealth from the Current Encounter menu and chooses a video display option.

Can display within Zedmed Clinical (ideal for one monitor) or in a web browser (ideal for two monitors).
6. The practitioner selects **Start Consult** in Telehealth.

The practitioner and patient see each other and can start talking.

## Considerations

The first time Telehealth is used on a new device or browser, you will be prompted to allow access to the camera and the microphone. Select **Allow** for this message. The request originates from the web browser used to facilitate the connection.

The patient should be **Admitted** from the Waiting Room to start the encounter, so the Admit time is recorded in the Waiting Room. If the patient is not Admitted, reception will not know if the patient is still waiting.

If the patient record is manually opened, a red banner will display at the top of the record when the patient joins the Telehealth session. Clicking the red banner will open the Telehealth session for the practitioner.

## How to start a Telehealth consult

To run the Telehealth consultation.


1. Monitor the **Waiting Room**.

Patients automatically attend to the Waiting Room when they open Telehealth on their phone and tap **Join**.

A **red** notification banner appears above the patient's record when they have selected **Join**.


2. Select **Admit**.

The patient's record opens in Zedmed and displays the **Current Encounter** modules, including **Telehealth**.



3. From the **Current Encounter** menu, open Telehealth using one of four options:


- **In Clinical** - opens Telehealth within Zedmed Clinical - **best option for one monitor**.
- **In Browsers** - opens Telehealth in your web browser - **best option for two monitors**.
- **Copy URL** - open your desktop browser and copy-paste the URL provided.
- **SMS Link to doctors** - sends a link to the doctor's phone using the mobile number in their Doctor Details.



4. Select **Start Consult**.

This displays the waiting patient and starts the Telehealth session.

The screenshot shows Telehealth running **In Clinical**.



If there is a disconnection, select **Re-Join** to resume the Telehealth consultation.

- When the leave Icon is selected, Telehealth will close. This does not end the Encounter in Zedmed.



## Patient workflow

Open the Telehealth link	Approve pictures and video	Tap Join
--------------------------	----------------------------	----------

