


# Appstream menu

Version: 1.00 | Last Modified on 27/10/2025 9:37 am AEDT

Some Zedmed customers use AppStream to access Zedmed Cloud because their practice setup or requirements differ from those using the standard interface. This guide explains the AppStream menu and its available functions.

## Appstream function menu



### A) Launch App

Use this menu if you need to restart Zedmed.

### (B) Switch windows - Switch between Clinical and Office

Select **Switch windows** to display Zedmed Office and Clinical side by side, then choose the one you want.

### C) My Files - upload to Zedmed


My Files is used to upload scans and documents to Zedmed Cloud using three intermediary folders.

- The **Zedmed Scans Folder** is a cloud folder that is linked to Clinical's Scans folder, which is used for incoming documents in Zedmed.
- The **Temporary Files** folder is used to upload a file from your PC to the cloud so you can then add the file to a patient's record.
- The **Home Folder** contains files that the user uploaded, and these files move with the user when they log in at different locations.

To learn more, see our dedicated Zedmed Cloud [File management guide](#).

### (D) Settings

Opens preferences for the streaming session. In Zedmed Cloud, the **Remember Audio and Video Settings** option determines whether the Microphone **mute** or **unmute** setting is remembered between sessions.



## (E) Microphone

Click on this icon to **mute** or **unmute** your microphone when using **Heidi AI** or **Telehealth**.

## (F) Fullscreen

Maximises the Zedmed Cloud window to full screen. You can hover your mouse at the top of the screen to view the menu bar. Press the **Esc** key on your keyboard to undo full-screen mode.

## (H) Notifications

Provides information about the streaming session, such as notifications about connections, network performance, or updates.

## (H) Profile - log out

To log out, select the Profile icon > **End session**.

To start a session again, select the Launch App as shown under item A in this guide.

## (I) Options

Additional settings and controls for the streaming session, such as the option to enter or exit full-screen mode, access the clipboard, or change display settings.