


Zedmed Mobile quickstart

Version: 1.00 | Last Modified on 26/11/2025 1:54 pm AEDT

Zedmed Mobile is **set up for the practitioner** and runs in your browser to view waiting patients and initiate Telehealth sessions. Before starting, log in to **Zedmed** then log in to **Zedmed Mobile** via your web browser:
<https://companion.tele-health.systems/>

1. A patient joins Telehealth and is auto-attended to the **Waiting Room** and shown in **Zedmed Mobile** as **green**.
2. When ready, **Admit** the patient from the **Waiting Room**, and they will appear in **Zedmed Mobile** as **orange**.



Step 1 - 'Admit' from Waiting Room
Step 2 - Click 'Open Telehealth'

Yellow: Patient joined + Dr Admitted
Green: Patient joined Telehealth
Grey: Patient has not joined


Waiting Room
Tuesday, 17 Jun 2025 **13:04 (-8) hours**

FILE #	Status	Patient Name	Appt Time	Attd Time	Admit Time	Doctor	Appt Type	Dept	Attd Notes
17	Admitted	Jones, Marrianna-Louise	13:00	12:51		PD	Telehealth		

Right-click an attendance to choose an action from the popup menu.

1 Admit

3. On Zedmed Mobile, select the **Open Telehealth** icon. This displays Telehealth on top of Zedmed Mobile.
4. Select **Start Consult**.



Practitioner's display

Attendees have arrived and are waiting

Welcome Dr Phillip Davis

Start Consult

Speaker Test

Patient's display

Marrianna-Louise Jones

Dr Phillip Davis

EXIT

