

# Telehealth Cloud consult

Version: | Last Modified on 04/12/2025 2:02 pm AEDT

Zedmed Cloud practitioners use Zedmed Mobile, an online portal, to facilitate Telehealth consultations. This portal runs in your web browser and allows you to easily monitor patients who have connected to Telehealth and are waiting, and to initiate Telehealth sessions.

Click this link for a **quick-start version of this guide** (a 1-page printable).

When you open the quickstart guide, select the PDF icon beside the quickstart page's title to download and print.

## Considerations

- The first time you use Zedmed Mobile, create a login and add your clinic as shown in the **Zedmed Mobile Setup guide**.
- When you use Zedmed, you should also open Zedmed Mobile, which authenticates via **Single Sign-On (SSO)**.
- Zedmed Mobile is opened in a web browser using the bookmarked URL: **<https://companion.telehealth.systems/>**
- When using a single monitor, position Zedmed Mobile in your browser at 1/4 or 1/3 screen width next to Zedmed's browser.

## Using Zedmed Mobile

Zedmed Mobile lets you see patients who have connected to Telehealth and are waiting for their consult - and start a telehealth consult with a waiting patient.

To perform a Telehealth consult with Zedmed Mobile:

1. Monitor the status of patients using the colour-codes below:
  - Patient has not tapped Join (practitioner cannot start Telehealth consult) = **Grey**.
  - Patient taps **Join**, which **Attends** them to the Waiting Room = **Green- Attended (Patient Arrived)**.
  - Practitioner **Admits** the patient from the Waiting Room = **Yellow - Admitted (Encounter Started)**.

Tuesday, 17 June PDAVIS 01:37:38 PM

[Show Billed](#) Hiding Billed Consults
[Telehealth](#) Filter
[Refresh](#)

Jones, Marrianna-Louise **Admitted** 01:00 PM **2** [Open Telehealth](#)  
Telehealth  
 Duration 15 mins

**Step 1 - 'Admit' from Waiting Room**  
**Step 2 - Click 'Open Telehealth'**

**Yellow:** Patient joined + Dr Admitted  
**Green:** Patient joined Telehealth  
**Grey:** Patient has not joined

Waiting Room

**Tuesday, 17 Jun 2025 13:04 (-8) hours**

17/06/2025 13:04 -8 Doctor Dr Phillip Davis Sort Attendance Time Filter Not Billed

FILE #	Status	Patient Name	Appt Time	Attd Time	Admit Time	Doctor	Appt Type	Dept	Attd Notes
17		Jones, Marrianna-Louise	13:00	12:51		PD	Telehealth		

< Previous Week < Previous Day Today > Next Day > Next Week

Right-click an attendance to choose an action from the popup menu.

[Refresh](#) **1** [Admit](#)

2. When ready, **Admit** the patient from the Zedmed **Waiting Room**.

This opens the patient's record in Clinical and starts the consultation timer.

3. On Zedmed Mobile, select the **Open Telehealth** icon (patient must be admitted-yellow)

Telehealth will open on the Zedmed Mobile webpage.

4. Select **Start Consult**.

Conduct the consultation as normal, using the video to view and communicate with the patient.

If there is a disconnection, select **Re-Join** to resume the Telehealth consultation.

5. Select **Exit** when the consultation has ended.

The patient will be removed from Zedmed Mobile once they have been billed.

← Back to Consults Telehealth 01:12:08 PM

**Practitioner's display**

Attendees have arrived and are waiting

Welcome Dr Phillip Davis

[Start Consult](#)

Speaker Test

< Test

← Back to Consults Telehealth 01:12:45 PM

**Patient's display**

Marrianna-Louise Jones

Dr Phillip Davis

[EXIT](#)

### Zedmed Mobile options:

- **Filters** - Choose between **All** appointments and **Telehealth** appointments only.
- **Tele Waiting** - Show only patients who have arrived (opened the link and clicked **Join**).

If they leave, they will change to **Telehealth** appointments.

- **Show Billed** - Includes patients who have been billed (The Waiting Room is checked for this status).
- **Refresh** - Zedmed Mobile refreshes the patients shown every minute, or when the **Refresh** button is selected.

## Switching Clinics

Zedmed Mobile supports multiple clinics, allowing practitioners to switch between each clinic's Waiting Room.

To change clinics:

1. Select the icon on the top right.
2. Select the Clinic from the list shown.

**Zedmed Mobile** will now show the **Waiting Room** for that Clinic.

To add another clinic, see the last section in the [Zedmed Mobile setup](#) guide.

