

Cannon Remote Capture


Version: 1.00 | Last Modified on 08/12/2025 11:31 am AEDT

Cannon Remote Capture allows Zedmed Cloud users to create custom folders where images can be saved, then imported into a patient's record. The images are imported in alphabetical order with the import time and date.


Setting up the capture folder

To create a practitioner-specifc folder:

- 1. In Clinical, select Tools > My Options.
- 2. In My options, locate Override image capture folder.
- 3. Tick Override.
- 4. Enter the folder name without spaces.
- 5. Select **OK** to save and exit.




After a few minutes, this folder will appear in My Files.



Using the capture folder

To save a patient's photos to their record using this folder:


1. Select My files, then upload or copy/paste the images into the folder you created during setup.



- 2. Open the patient in **Clinical** and start the encounter.
- 3. Select Capture > Canon Remote Capture.



4. The **Canon Remote Capture** folder will open, containing the files uploaded in step 1.



- 5. Select **Import All**, and the images will be uploaded to the patient's record.
- 6. Open Summary Views > Images/ECG to view the images.

